





# Acetyl-UBA52/RPS27A/UBB/UBC (K29) Antibody

<b>Product Code</b>	CSB-PA000143
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P62987/P62979/P0CG47/P0CG48
<b>Immunogen</b>	Synthesized peptide derived from the Human Ub around the acetylation site of K29.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	UBA52/RPS27A/UBB/UBC

## Image



Western Blot analysis of HepG2 cells using Acetyl-Ub (K29) Polyclonal Antibody



Western Blot analysis of HepG2 cells using Acetyl-Ub (K29) Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.