




Acetyl-HMGB1 (K12) Antibody

Product Code	CSB-PA000149
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09429
Immunogen	Synthesized peptide derived from the N-terminal region of Human HMG-1 around the acetylation site of K12.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HMGB1; HMG1; High mobility group protein B1; High mobility group protein 1; HMG-1
Immunogen Species	Homo sapiens (Human)
Target Names	HMGB1

Image



Western Blot analysis of HEPG2-UV cells using Acetyl-HMG-1 (K12) Polyclonal Antibody