




Strep-Tag Monoclonal Antibody

Product Code	CSB-MA000171
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	N/A
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Tag Control Antibody
Immunogen Species	N/A
Clone No.	19B4

Image



1µg Strep II fusion protein+ Primary antibody dilution at 1) 1:50002) 1:10000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.