




# Myc-Tag Monoclonal Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-MA000179   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Immunogen</b>           | Synthetic Peptide  |
| <b>Raised In</b>           | Mouse  |
| <b>Species Reactivity</b>  | N/A  |
| <b>Tested Applications</b> | ELISA,WB;Recommended dilution:WB:1:500-1:5000                                |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | 0  |
| <b>Product Type</b>        | Monoclonal Antibody  |
| <b>Species</b>             | N/A  |

## Image



Western blot analysis of Recombinant protein, diluted at 1) 1:5000 2) 1:10000