




# GAPDH Monoclonal Antibody

<b>Product Code</b>	CSB-MA000184
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04406
<b>Immunogen</b>	Recombinant full-length protein
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	H,R,M,Mk,Dg,C,Hm,Rb,Pg,Sh,Insect,Yeast
<b>Tested Applications</b>	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
<b>Isotype</b>	IgG
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	H, R, M, Mk, Dg, C, Hm, Rb, Pg, Sh, Insect, Yeast
<b>Target Names</b>	GAPDH
<b>Clone No.</b>	11H4

## Image



Western blot analysis of Hela (1), Rat brain (2), Rabbit Muscle (3), Sheep Muscle (4), and Mouse brain (5), diluted at 1:10000.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.