





PCNA Monoclonal Antibody

Product Code	CSB-MA000186
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12004
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% BSA.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	H, R, M
Target Names	PCNA
Clone No.	19F7

Image



IHC staining of Human breast cancer tissue paraffin-embedded, diluted at 1:200



Western blot analysis of HeLa (1), Rat brain (2), NIH 3T3 (3), 293T (4), diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.