





PCNA Monoclonal Antibody

Product Code	CSB-MA000186
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12004
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	PCNA; Proliferating cell nuclear antigen; PCNA; Cyclin
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	PCNA

Image



IHC staining of Human breast cancer tissue paraffin-embedded, diluted at 1:200



Western blot analysis of Hela (1), Rat brain (2), NIH 3T3 (3), 293T (4), diluted at 1:5000.