




ACTB Monoclonal Antibody

Product Code	CSB-MA000187
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	H,R,M,Mk,Dg,Ch,Hm,Rb,Pg,Sh
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	ACTB; Actin, cytoplasmic 1; Beta-actin
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB
Clone No.	1G8

Image



Western blot analysis of HeLa (1), Rat brain (2), Mouse brain (3), (A02) diluted at 1:10000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.