




ACTB Monoclonal Antibody

Product Code	CSB-MA000187
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	H,R,M,Mk,Dg,Ch,Hm,Rb,Pg,Sh
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	ACTB; Actin, cytoplasmic 1; Beta-actin
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	ACTB

Image



Western blot analysis of HeLa (1), Rat brain (2), Mouse brain (3), (A02) diluted at 1:10000.