




COX4I1 Monoclonal Antibody

Product Code	CSB-MA000188
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13073
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	COX4I1; COX4; Cytochrome c oxidase subunit 4 isoform 1, mitochondrial; Cytochrome c oxidase polypeptide IV; Cytochrome c oxidase subunit IV isoform 1; COX IV-1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COX4I1
Clone No.	30H2

Image



Western blot analysis of HeLa, diluted at 1) 1:2000 2) 1:5000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.