





TUBA1A Monoclonal Antibody

Product Code	CSB-MA000190
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q71U36/P68363
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IP;Recommended dilution:WB:1:500-1:5000,IP:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	Tubulin alpha-1A chain
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	10C8

Image



1) Input: Mouse Brain Tissue Lysate 2) IP product: IP dilute 1:200



Western blot analysis of 1) HeLa, 2) Rat Brain Tissue, 3) Mouse Brain Tissue, diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.