




GAPDH Monoclonal Antibody

Product Code	CSB-MA000195
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04406
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	H,R,M,Mk,Dg,C,Hm,Rb,Pg,Sh
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GAPDH
Clone No.	13B6

Image



Western blot analysis of Hela, diluted at 1:3000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.