





HSPA1L Monoclonal Antibody

Product Code	CSB-MA000197
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34931/P08107
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Heat shock 70 kDa protein 1-like/Heat shock 70 kDa protein 1A/1B
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	21E12

Image



IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.



Western blot analysis of 1) HeLa, 2) Mouse Brain, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.