





LGALS3 Monoclonal Antibody

Product Code	CSB-MA000198
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17931
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB: 1:500-1:5000, IHC: 1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	LGALS3; MAC2; Galectin-3; Gal-3; 35 kDa lectin; Carbohydrate-binding protein 35; CBP 35; Galactose-specific lectin 3; Galactoside-binding protein; GALBP; IgE-binding protein; L-31; Laminin-binding protein; Lectin L-29; Mac-2 antigen
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LGALS3
Clone No.	11F9

Image



IHC staining of Human thyroid tissue paraffin-embedded, diluted at 1:200.



Western blot analysis of Hela, diluted at 1:3000.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.