




PTPRC Monoclonal Antibody

Product Code	CSB-MA000202
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08575
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB: 1:500-1:5000, IHC 1:50-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	PTPRC; CD45; Receptor-type tyrosine-protein phosphatase C; Leukocyte common antigen; L-CA; T200; CD45
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTPRC
Clone No.	5B12

Image



Western blot analysis of Hela, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.