




# AMACR Monoclonal Antibody

<b>Product Code</b>	CSB-MA000203
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UHK6
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	AMACR; Alpha-methylacyl-CoA racemase; 2-methylacyl-CoA racemase
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Gene Names</b>	AMACR

## Image



IHC staining of mouse prostate adenocarcinoma tissue, diluted at 1:200.



Western blot analysis of 1) HepG2, 2) Mouse Kidney, diluted at 1:1000.