




AMACR Monoclonal Antibody

Product Code	CSB-MA000203
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHK6
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	H, M, R
Target Names	AMACR
Clone No.	7A7

Image



IHC staining of mouse prostate adenocarcinoma tissue, diluted at 1:200.



Western blot analysis of 1) HepG2, 2) Mouse Kidney, diluted at 1:1000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.