




FCER2 Monoclonal Antibody

Product Code	CSB-MA000205
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06734
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	FCER2; CD23A; CLEC4J; FCE2; IGEBF; Low affinity immunoglobulin epsilon Fc receptor; BLAST-2; C-type lectin domain family 4 member J; Fc-epsilon-RII; Immunoglobulin E-binding factor; Lymphocyte IgE receptor; CD23
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCER2
Clone No.	12C2

Image


IHC staining of Human tonsil tissue paraffin-embedded, diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.