




# GFAP Monoclonal Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-MA000216   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Uniprot No.</b>         | P14136   |
| <b>Immunogen</b>           | Synthetic Peptide  |
| <b>Raised In</b>           | Mouse  |
| <b>Species Reactivity</b>  | Rat,Mouse  |
| <b>Tested Applications</b> | ELISA,WB;Recommended dilution:WB:1:500-1:5000                                |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | GFAP; Glial fibrillary acidic protein; GFAP                                  |
| <b>Product Type</b>        | Monoclonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Gene Names</b>          | GFAP   |

## Image



Western blot analysis of Rat Brain Tissue, diluted at 1:5000.