




GFAP Monoclonal Antibody

Product Code	CSB-MA000216
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14136
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Rat,Mouse
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	GFAP; Glial fibrillary acidic protein; GFAP
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GFAP
Clone No.	13E9

Image



Western blot analysis of Rat Brain Tissue, diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.