




CEACAM5 Monoclonal Antibody

Product Code	CSB-MA000218
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06731
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	CEACAM5; CEA; Carcinoembryonic antigen-related cell adhesion molecule 5; Carcinoembryonic antigen; CEA; Meconium antigen 100; CD66e
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEACAM5
Clone No.	6G12

Image



IHC staining of human colon cancer tissue, diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.