




CEACAM5 Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA000218 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P06731 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human |
| Tested Applications | ELISA,IHC;Recommended dilution:IHC:1:100-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | CEACAM5; CEA; Carcinoembryonic antigen-related cell adhesion molecule 5; Carcinoembryonic antigen; CEA; Meconium antigen 100; CD66e |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Gene Names | CEACAM5 |

Image



IHC staining of human colon cancer tissue, diluted at 1:200.