




MKI67 Monoclonal Antibody

Product Code	CSB-MA000222
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46013
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	H, M, R
Target Names	MKI67
Clone No.	19A10

Image


IHC staining of Mouse Kidney tissue, diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.