





CTNNB1 Monoclonal Antibody

Product Code	CSB-MA000223
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	CTNNB1

Image



IHC staining of human breast cancer tissue, diluted at 1:200.



Western blot analysis of 1) Hela, 2) 293T, 3) MCF7, 4) Mouse Brain tissue, 5) Rat LiverTissue, diluted at 1:2000.