



CD2 Monoclonal Antibody

Product Code	CSB-MA000227
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06729
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC; Recommended dilution: IHC: 1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	CD2; SRBC; T-cell surface antigen CD2; Erythrocyte receptor; LFA-2; LFA-3 receptor; Rosette receptor; T-cell surface antigen T11/Leu-5; CD2
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD2
Clone No.	20B8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.