



MME Monoclonal Antibody

Product Code	CSB-MA000230
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08473
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	MME; EPN; Nephilysin; Atriopeptidase; Common acute lymphocytic leukemia antigen; CALLA; Enkephalinase; Neutral endopeptidase 24.11; NEP; Neutral endopeptidase; Skin fibroblast elastase; SFE; CD10
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MME
Clone No.	24E3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.