




FH Monoclonal Antibody

Product Code	CSB-MA000231
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07954
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Fumarate hydratase, mitochondrial
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	FH

Image



Western blot analysis of 1) 293T, 2) HepG2, 3) Hela, diluted at 1:3000.