




ABCB5 Monoclonal Antibody

Product Code	CSB-MA000232
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2M3G0
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	ABCB5; ATP-binding cassette sub-family B member 5; ABCB5 P-gp; P-glycoprotein ABCB5
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCB5
Clone No.	12F4

Image



Western blot analysis of 1) HeLa, 2) 293T, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.