




FCGR3A Monoclonal Antibody

Product Code	CSB-MA000248
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08637/O75015
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Low affinity immunoglobulin gamma Fc region receptor III-A/Low affinity immunoglobulin gamma Fc region receptor III-B
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	22C4

Image



Western blot analysis of 1) Jurkat, 2) K562, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.