




# Mono-Methyl-Histone H3(K79) Monoclonal Antibody

|                            |                                                                                                                                                                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-MA000251                                                                                                                                                                                                                                                    |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                   |
| <b>Uniprot No.</b>         | P68431/Q71DI3/P84243                                                                                                                                                                                                                                            |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                                                                                                               |
| <b>Raised In</b>           | Mouse                                                                                                                                                                                                                                                           |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                                                                                                                                                                 |
| <b>Tested Applications</b> | ELISA,WB;Recommended dilution:WB:1:500-1:5000                                                                                                                                                                                                                   |
| <b>Form</b>                | Liquid                                                                                                                                                                                                                                                          |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                                                                                                                                  |
| <b>Storage Buffer</b>      | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.                                                                                                                                                                           |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                             |
| <b>Alias</b>               | HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/ |
| <b>Product Type</b>        | Monoclonal Antibody                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                            |
| <b>Clone No.</b>           | 3E10                                                                                                                                                                                                                                                            |

## Image



Western blot analysis of 1) Hela, 2) Rat Testis tissue, 3) Raw264.7, diluted at 1:2000.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.