




KRT18 Monoclonal Antibody

Product Code	CSB-MA000255
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05783
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	KRT18; CYK18; PIG46; Keratin, type I cytoskeletal 18; Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	KRT18

Image



Western blot analysis of 1) Hela, 2) HepG2, 3) Mouse Skeletal Muscle tissue, (2F7) diluted at 1:1000.