



FUT4 Monoclonal Antibody

Product Code	CSB-MA000264
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22083
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	FUT4; ELFT; FCT3A; Alpha-(1,3)-fucosyltransferase; ELAM-1 ligand fucosyltransferase; Fucosyltransferase 4; Fucosyltransferase IV; Fuc-TIV; FucT-IV; Galactoside 3-L-fucosyltransferase
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FUT4
Clone No.	25G5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.