



ITGA2B Monoclonal Antibody

Product Code	CSB-MA000265
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08514
Immunogen	Synthetic Peptide of CD41
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	Human
Gene Names	ITGA2B