





RELA Monoclonal Antibody

Product Code	CSB-MA000270
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04206
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IP;Recommended dilution:WB:1:500-1:5000,IP:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	RELA; NFKB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RELA
Clone No.	26E8

Image



1) Input: HeLa Cell Lysate 2) IP product: IP dilute 1:200



Western blot analysis of 1) HeLa, 2) Rat Heart Tissue, 3) Mouse Spleen Tissue, diluted at 1:2000.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.