





PRDX1 Monoclonal Antibody

Product Code	CSB-MA000271
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06830
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	PRDX1; PAGA; PAGB; TDPX2; Peroxiredoxin-1; Natural killer cell-enhancing factor A; NKEF-A; Proliferation-associated gene protein; PAG; Thioredoxin peroxidase 2; Thioredoxin-dependent peroxide reductase 2
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRDX1
Clone No.	24F11

Image



IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.



Western blot analysis of 1) MCF7, 2) Rat Kidney Tissue, 3) Mouse Brain Tissue, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.