




TUBB3 Monoclonal Antibody

Product Code	CSB-MA000274
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13509
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUBB3
Clone No.	4D4

Image



Western blot analysis of Zebrafish skeletal muscle, (Zebrafish Specific) diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.