



Insulin Glargine Monoclonal Antibody

Product Code	CSB-MA000275
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Insulin Glargine
Raised In	Mouse
Species Reactivity	All
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	N/A
Clone No.	13E8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.