




Di-methyl-Histone H3(K79) Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA000276 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P68431/Q71DI3/P84243 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/ |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Clone No. | 3F8 |

Image



Western blot analysis of 1) HeLa, 2) Raw264.7, 3) Rat Brain Tissue, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.