




TUBB1 Monoclonal Antibody

Product Code	CSB-MA000278
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4B7
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	TUBB1; Tubulin beta-1 chain; TUBB2A; TUBB2; Tubulin beta-2A chain; TUBB2B; Tubulin beta-2B chain; TUBB4B; TUBB2C; Tubulin beta-4B chain; Tubulin beta-2 chain; Tubulin beta-2C chain; TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-II
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUBB1
Clone No.	8H4

Image



Western blot analysis of Zebrafish skeletal muscle, diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.