




ACTB Monoclonal Antibody

Product Code	CSB-MA000280
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% BSA
Isotype	IgG
Alias	ACTB; Actin, cytoplasmic 1; Beta-actin
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB
Clone No.	9A9

Image



Western blot analysis of Zebrafish skeletal muscle, diluted at 1:3000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.