




COL3A1 Monoclonal Antibody

Product Code	CSB-MA000282
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02461
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	COL3A1; Collagen alpha-1(III) chain
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL3A1
Clone No.	16C2

Image



Western blot analysis of Hela, diluted at 1:1000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.