





GFP Monoclonal Antibody

Product Code	CSB-MA000283
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	All
Tested Applications	ELISA,WB,IP;Recommended dilution:WB:1:500-1:5000,IP:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	N/A
Clone No.	9D11

Image



1) Input (control) 2) IP products, antibody dilution 1:200



Western blot analysis of GFP transfected HeLa, diluted at 1) 1:5000 2) 1:10000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.