




PARP1 Monoclonal Antibody

Product Code	CSB-MA000290
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09874
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Poly [ADP-ribose] polymerase 1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PARP1
Clone No.	10G12

Image



Western blot analysis of 1) HeLa, 2) 293T, 3) Jurkat, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.