




ACTB Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA000292 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P60709 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human,Mouse,Rat,Mk,Dg |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA |
| Isotype | IgG |
| Alias | ACTB; Actin, cytoplasmic 1; Beta-actin |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ACTB |
| Clone No. | 6B6 |

Image



Western blot analysis of 1) HeLa, 2) Mouse Brain tissue, 3) Rat Brain tissue, diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.