





EGFR Monoclonal Antibody

Product Code	CSB-MA000293
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00533
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EGFR
Clone No.	24C3

Image



Western blot analysis of 1) A431, 2) HeLa, diluted at 1:2000.



1?Input: HeLa Cell Lysate 2?IP product: IP dilute 1:200 Western blot analysis: primary antibody : 1:1,000 Secondary antibody: Goat anti-Mouse IgG, Heavy Chain Specific(S003), 1:5,000



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.