





EIF4A1 Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA000295 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P60842 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | Eukaryotic initiation factor 4A-I |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Gene Names | EIF4A1 |

Image



IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.



Western blot analysis of 1) 293T, 2) HeLa, 3) HepG2, 4) Mouse Brain tissue,