





# EIF4A1 Monoclonal Antibody

<b>Product Code</b>	CSB-MA000295
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P60842
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	H, M, R
<b>Target Names</b>	EIF4A1
<b>Clone No.</b>	30E9

## Image



IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.



Western blot analysis of 1) 293T, 2) HeLa, 3) HepG2, 4) Mouse Brain tissue,

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.