




FN1 Monoclonal Antibody

Product Code	CSB-MA000296
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02751
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FN1
Clone No.	29E4

Image



Western blot analysis of HeLa, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.