




His-Tag Monoclonal Antibody

Product Code	CSB-MA000300
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	N/A
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0
Product Type	Monoclonal Antibody
Immunogen Species	N/A
Target Names	His-Tag
Clone No.	8B3

Image



Western blot analysis of His-Tag recombinant protein, diluted at 1:10000 1) His-MBP recombinant protein 2) His-Tag transfected 293T

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.