




PARP1 Monoclonal Antibody

Product Code	CSB-MA000305
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09874
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	PARP1; ADPRT; PPOL; Poly [ADP-ribose] polymerase 1; PARP-1; ADP-riboseyltransferase diphtheria toxin-like 1; ARTD1; NAD(+) ADP-riboseyltransferase 1; ADPRT 1; Poly[ADP-ribose] synthase 1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PARP1
Clone No.	14B6

Image



Western blot analysis of Jurkat, diluted at 1:3000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.